

Stannah

Stratum Homelift User Manual

ENG Original instructions



CONTENTS

Description	4
General Dos and Don'ts	5
Controls and Operation	7
Smoke Alarms	9
Emergency Procedures	11
Fault Finding	14
Changing Call Station Batteries	16
Lift Disassembly / Safe Disposal of Hazardous Materials	18
Service History Record	19
Declaration of Conformity	20
Lift Specification	21

DESCRIPTION

The Stratum Homelift is an inter floor lift that is designed for use by persons with impaired mobility travelling between fixed floor levels in private dwellings.

With a maximum carrying capacity of one person, with or without a wheelchair (or 250kg), this lift is not intended for use as a means of transporting goods.

The lift is designed to operate without a lift shaft and is provided with an automatic in fill panel which makes the ceiling aperture safe when the lift is parked downstairs.

Two way communication system should be specified for emergency communication.

A standard feature is the provision of half hour fire rated panels in both the aperture in fill and the car underpan.

The lift car panels are made from powder coated steel which can easily be cleaned using normal household cleaners. Upholstery is made from PVC and can be cleaned in the same way.

A smoke detection system has been installed on your lift. It has been designed to provide adherence to British Standard BS5900 2012 Section 9.13 “Behaviour of homelift in the event of fire”.

GENERAL DOS AND DON'TS

- Never switch off the power supply to the lift, even when you go away. The lift control circuits are fed by a battery, which must be kept on constant charge.
- The lift should always be returned to the lower level when not in use. If it is left upstairs for prolonged periods, it will occasionally re-level itself depending on conditions. The lift must be left at the lower level if the mains is turned off.
- If your lift is fitted with a manual door always close it after use. Powered door units have a self closer. Do not pull or push the automatic door.
- Never allow children to play in, under or around the lift. If children are in the house, isolate the lift using the optional remote isolate fob, (page 7).
- Ensure that the area under the lift is kept clear. The underpan safe surface is fitted with sensors, which automatically stop the lift if it strikes an object (page 6).
- Always keep your emergency door release key and optional remote isolate fob in the lift or in a safe place near the lift.
- Do not place any object on the aperture infill panel safe edge or stand on it when the lift is in operation. Ensure that as far as practical, the area around the travelling infill panel safe edge is clear of persons (particularly children) when the lift is being operated. This is to ensure there is no danger of them falling into the car when the lift is in use. The infill panel safe edge is fitted with sensors that automatically stop the lift if the infill panel is obstructed (page 6)
- Do not use this lift for anything other than transporting those with impaired mobility between fixed floor levels
- Do not lean over the car sides or door. These are fitted with safety edges which will stop the lift if activated.
- Always treat your lift with respect that should be shown to electrical and mechanical equipment.

- Users with wheelchairs should apply the brakes on their chair and all other users should ensure they are using the seat provided before moving the lift. Do not travel in the lift unless seated.
- Safety related components should only be adjusted and reset by a competent person.
- Proceed with caution when exiting the carriage backwards.

The diagram below (*Figure 1*) shows the position of all the sensors on the lift, designed to prevent injury or damage if the movement of the carriage is obstructed.




Figure 1. lift shown part-way through aperture

CONTROLS AND OPERATION

There are two wall mounted call stations, one at each level (Figure 2) and a similar control station fitted in the lift car. The lift can be isolated by using the optional remote isolate fob (Figure 3). When the lift is isolated, none of the control stations will function. The call and control stations can only be activated by using the optional remote isolate fob. When the lift is activated, the coloured indicator lights in the car will illuminate. The lights in the car will switch on automatically when any call or control button is pressed and will automatically turn off after a few minutes.

Once the lift has been stopped it cannot be restarted for 3 seconds.




Figure 2. Wall mounted call station

Figure 3. Optional remote isolate fob

Users with wheelchairs should apply the brakes on their chair and all other users should ensure they are using the seat provided before starting the lift.

General operation

Call the lift by pressing the up arrow or down arrow button on either call station and wait for it to stop. The door will automatically open. Enter the lift and press and release the door button to close it. Then press and release UP or DOWN on the control station. The lift will travel uninterrupted to the next floor. If the lift does not start, check that the door is properly closed and try again. Always close the door after using the lift by pressing the door button on the wall mounted call station. Leave the lift downstairs whenever possible.

In the event that the normal controls fail to work, emergency controls can be found under the flap above the mirror. Green is 'UP', red is 'STOP', black is 'DOWN' and white is 'DOOR' (Figure 4).




Figure 4. Wall mounted call station




Figure 5. Optional remote isolate fob

When outside the car, in an emergency the door can be opened by pressing the front end of the side safe edge, on the corner where the door opens (Figure 5).

SMOKE ALARM

A smoke detection system has been installed on your lift. It has been designed to provide adherence to British Standard BS5900 2012 Section 9.13 “Behaviour of homelift in the event of fire”. The system utilises two smoke alarms one upper level, one lower level, which are wirelessly connected to the main circuit board on the lift (Figure 6).




Figure 6.

The Smoke Alarm unit comprises two parts, the first (TOP PART) being a Smoke alarm unit. This is the Smoke detection unit and contains an integral battery with a 10 year life span. This is not replaceable.

The second part (BASE UNIT) contains the radio transmitter receiver unit. Power to this unit is supplied by single 10 year life span lithium battery.

Radio smoke alarm operation

When installed on a Harmony Lift, the smoke alarm system will cause the lift to deactivate safely once the alarm is triggered. When the lift is deactivated, the door will continue to operate as normal.

When the Lift is stationary at either level

If smoke is detected, the alarm will sound. After a period of time, all other smoke alarms connected to the system will then start to sound and the lift deactivates.

When the Lift is travelling between levels

If smoke is detected, the alarm will sound. After a period of time, all other smoke alarms connected to the system will then start to sound.

The lift will continue to its requested level, it will remain possible (until that level is reached) to change the direction of the lift. Once at the desired level, the lift will deactivate.

Reactivation of Lift

The lift will automatically reactivate when the smoke alarm no longer detects smoke and a period of two minutes has expired.

Silencing the Smoke Alarms

The alarms can be silenced by pressing the mute button on the sensor that initiated the alarm. The initiating sensor can be identified by a Red light flashing every second.

When the alarms are deactivated, the lift will automatically reactivate when a period of two minutes has expired. If the source of the smoke alarm is not removed, the smoke alarms will begin to sound again and the lift will be disabled.

EMERGENCY PROCEDURES

In the event of a mains failure during travel, the battery backed control system of the lift will only allow normal operation in the down direction without loss of any safety features. This allows the user to exit the car at the lower level in the normal way.

Emergency Manual Lowering

IMPORTANT:

During emergency manual lowering, the normal safety features will not function, so the lift will not stop if a person, pet or object is under the lift.

- The exact lowering procedure must be observed, because the normal safety features will not function during manual lowering.
- The emergency lowering procedures should never be used if the lift is fully up or no one is trapped in it.
- The emergency lowering procedures should never be used as the normal down travel function until an engineer attends.

PLEASE NOTE:

If the lift is fully upstairs and a person is trapped inside, please see the section Emergency unlocking on pages 12-13

If the lift has stopped mid level and the customer is unable to get the lift up or down then the only time it should be lowered by the emergency valve is if:

- There is a second person (B) around the lift area at the lower level to ensure that nothing goes under the lift during the lowering by the first person (A)
- OR the person lowering the lift has sight of the area under the lift

Person A

- Ensure the door is fully closed
- Turn the mains supply to the lift off.
- Locate the hydraulic power unit, normally outside the property (Figure 7)

- Swivel the metal cover plate on the front face of the housing to reveal an access hole (Figure 7)
- The red cord revealed in the access hole now needs to be pulled continuously to lower the lift car slowly.
- After 5 seconds release cord and check with person B that the aperture infill panel is following the carriage. If so, resume pulling the cord as before.
- Once the lift is at the lower level turn the mains supply back on and secure metal cover plate.




Figure 7.

Person B

Remains in the house by the lift and communicates with Person A to ensure the safe lowering of the lift.

- Ensure that no object, person or pet are in the path of the lift travel
- Confirms that the trap door follows the lift during descent and locates fully in the floor to guard against the possibility of anyone falling down the lift way.

Emergency unlocking

WARNING: RISK OF FALLING – Emergency unlocking must only be undertaken when the lift is at the upper or lower landing level.

The lift car door is designed so that it will only open when the lift is at or within 25mm of each floor served. If for any reason the door cannot be opened, the door lock can be over-ridden from the inside using the following steps:

- a Remove the black rubber grommet by the door
- b Insert the black emergency door release key into the square hole behind the grommet
- c Turn the key to release the mechanism

This should release the catch. The door can be opened from the outside in the same way.

The electrically operated door on the wheelchair model can be forced open or closed manually. This will disengage the door from the opening mechanism and require re-setting of the door. Re-setting of the door can be done as follows:

- a Manually close the door as far as possible
- b Press and release the door button
- c Move the door back and forth until it clicks
- d Press and release the door button again

The door should now operate normally, otherwise repeat the above process.

FAULT FINDING

The most likely causes of your lift failing to operate are:

- The door not being fully closed
- Something obstructing the travel of the car or infill panel safe edge causing one or more of the sensors or safe edges (Figure 1) to operate and therefore preventing travel in the appropriate direction.

To assist in identifying the cause the car is fitted with a simple system of coloured indicator lights on the rear panel, one red, two yellow and two green, (Figure 8).




Figure 8.

Lift Malfunction

Fault	Indication	Cause	Remedy
Lift will not travel in either direction	No lights on car panel	Door not shut or remote fob off	Press door button Press a button on fob
		Flat batteries in smoke alarms	Replace smoke alarm batteries
Lift will not go up	One green off on car panel	Car safe edge obstruction	Remove obstruction or free safe edge
	Both greens off on car panel	Infill panel obstructed on upper surface	Remove obstruction from upper surface
Lift will not go down	One yellow off on car panel	Car underpan obstruction	Remove obstruction from beneath surface
	Both yellow off on car panel	In fill panel obstructed on lower surface	Remove obstruction from lower surface
Powered door will not close fully	Can be moved easily by hand	Door has disengaged from drive mechanism	See page 12-13 Emergency unlocking

Smoke Alarm Malfunction

Indication	Cause	Remedy
Alarm beeps every 11 minutes	Test Button Jammed	Remove from ceiling and release button
Alarm beeps twice every minute	Unit Malfunction	Call an engineer

Low Battery Warning

Indication	Cause	Remedy
Single short beep on lift car every 2 minutes	Low Handset battery	Replace CR2450 batteries in handset
Smoke Alarm beeps every minute	Low Smoke Alarm Battery	Replace unit

CHANGING CALL STATION BATTERIES




Figure 9.




Figure 10.




Figure 11.

1. The call station can be removed from the wall by sliding the case upwards.
2. Remove the four screws in the back with a posi-drive screwdriver.
3. Using a screw driver, gently push the battery out a short way and then pull using fingers.

NOTE: It is critical that all three batteries are replaced with new ones of the same type, manufacture and age, that they are fitted at the same time and that they are correctly oriented.

4. Refit back cover and re-attach to wall.

NOTE: If recessed controls are fitted, remove module by undoing two wall fixing screws then follow instructions above. Care must be taken not to damage protruding aerial during this process.

LIFT DISASSEMBLY/SAFE DISPOSAL OF HAZARDOUS MATERIALS

This lift must be disassembled by a competent person who has been fully trained in the installation of this lift and is qualified to provide safe disconnection of the lift to the mains terminal.

Batteries & Printed Circuit Boards (PCB)

The batteries and PCB's within this product should not be disposed of with other household waste at the end of their working life. Where batteries are marked with the chemical symbols Hg, Cd or Pb, it indicates that the battery contains mercury, cadmium or lead above the reference levels in EC Directive 2006/66. If batteries are not properly disposed of, these substances can cause harm to human health or the environment.

Batteries and PCB's that are no longer required for this lift, at the end of their working life, can be returned either to an approved waste disposal facility or to Terry Group Ltd for safe disposal.

Oil

Oil from this lift should be disposed of via an authorised waste disposal contractor, to an approved waste disposal facility.

SERVICE HISTORY RECORD

An entry should be added to the following table every time the lift is serviced

Date	Engineer	Company	Comments

DECLARATION OF CONFORMITY



Machinery Description: Stratum Homelift

This lift was manufactured by TERRY GROUP Ltd., who declare that this lift fulfils all the relevant provisions of the following Directives:

2014/30/EU	Electromagnetic Compatibility Directive
2006/42/EC	Machinery Directive

This lift also fulfils all the relevant provisions of the following Standards:

BSEN 12015:2014	Electromagnetic compatibility. Product family standard for lifts, escalators and moving walks. Emission.
BSEN 12016:2013	Electromagnetic compatibility. Product family standard for lifts, escalators and moving walks. Immunity.
BS5900:2012	Powered homelifts with partially enclosed carriers and no liftway enclosures – Specification

This Declaration of Conformity covers all lifts with serial numbers starting with L, H and R.

Person authorised to compile Technical File:

Greg Gnyp, Terry Group Ltd., Longridge Trading Est, Knutsford, Cheshire, WA16 8PR

EC examination carried out by: Bureau Veritas UK Ltd., Parklands, Wilmslow Road, Didsbury, Manchester, M20 2RE.

Notified Body Reference Number:0041

EC examination certificate number: CE-0041-MD-TER004-18-GBR

This declaration was completed at Terry Group Ltd., Longridge Trading Estate, Knutsford, Cheshire, WA16 8PR, in April 2018.

This compliance is only valid if the installation test Certificate has been completed and signed by a competent lift engineer trained to install this product to the latest installation instructions.

TERRY GROUP Ltd.

P.Morrey (Managing Director)

LIFT SPECIFICATION

Address of manufacturer:-

Terry Group Ltd.

Unit 1 Longridge Trading

Estate

Knutsford

Cheshire

England

WA16 8PR

Lift serial No: Year of manufacture:

Safe working load 250kg

Maximum travel 3.6 metres

Duty cycle 10 cycles per hr with max load

Average noise level 65 dBA

Power supply Dedicated 240V ~ 50/60 Hz
single phase supply

Control voltage 12V DC

Hydraulic pump power
consumption 750W maximum

Hydraulic oil grade T22

Hydraulic pump enclosure IP54

Test specification Manufactured and tested to
BS5900:2012

Fire specification Load-bearing capacity
30 minutes, integrity and
insulation 15 minutes.
Assessed by Warrington Fire
Research Centre.

For technical help,
sales or service enquiry telephone:



Stannah

Stannah Homelifts
Watt Close, East Portway,
Andover, Hampshire, SP10 3SD

t: 0800 715 482
e: webenquires@stannah.co.uk
www.stannahhomelifts.co.uk